

# MTS / T-BERD 6000A V2

# B15312

# RELEASE NOTE

## **Ste-Core Package 15312 (July 2015)**

- Boot V5.84 (14/July/2015)
- Linux Kernel V7.14 (28/July/2015)
- File System V5.32 (15/July/2015)
- Instrument Setup V14.50 (17/July/2015)
- Microscope V14.50 (17/July/2015)

## **Ste-Optical Package 15312 (July 2015)**

- Fiber Optics V14.50 (17/July/2015)

# CONTENTS

T-BERD/MTS 6000AV2 Base Unit.....	3
Improvements .....	3
OTDR .....	3
New functionalities .....	<b>Erreur ! Signet non défini.</b>
Improvements .....	3
OSA.....	3
Improvements .....	3
HR-OSA .....	3
Improvements .....	3
ODM.....	3
Improvements .....	3
Scope.....	3
New functionalities .....	3
Improvements .....	3
Script.....	4
Improvements .....	4

## T-BERD/MTS 6000AV2 Base Unit

### *Improvements*

- Improved pointer management
- Improved management of capacitive screen (less sensitivity)
- Improved management of screen saver
- Multitouch zoom deactivated with capacitive screen
- Improved management of pass/fail with Stratasync function

## OTDR

### *Improvements*

- Improved management of MM/SM module when moving from one tab to another
- Better management of configuration files for FTTA Extended-SLM mode

## OSA

### *Improvements*

- Improved management of OSA acquisition having deep rising and descending spectral shapes (was generating ISU-323 error)

## HR-OSA

### *Improvements*

- The minimum channel spacing can now be down to 0.5GHz

## ODM

### *Improvements*

- Improved tuning of ODM module

## Scope

### *New functionalities*

- Addition of comments into scope report

### *Improvements*

- Improved auto increment during storage of image + PDF files
- Improved auto-save with P5000 scope

- Improved reading of scope image readings from the explorer of the scope function
- The scope button is reactivated and acts as “Test” softkey on platform

---

## Script

---

### *Improvements*

- Improved keyboard management
- Improved close of script when switching off the platform while FCHAR function is active
- The fiber information fields is enabled with up to 40 characters
- Improved management of the script configuration file